

EduPrimer Alu DNA Profiling Kit

Kit Components

<u>Kit Components Included:</u>	<u>Qty</u>	
DNA Isolation Components		
--Solution A (Lysis buffer)	6 mL	included
--Solution B (Neutralizer)	0.6 mL	included
--Cotton Swabs	28 swabs	included
PCR Reaction Components		
--2X PCR Master Mix (includes TPA25 or PV92 primers)	300 μ L	included
--Positive Control DNA(Heterozygous)	30 μ L	included
-- H ₂ O/ Negative Control	300 μ L	included
Gel Electrophoresis Components		
--DNA Ladder	30 μ L	included