

# Animal Identification: Catch a Predator Kit I and II

## Kit Components

| <b><u>Kit Components Included:</u></b>                                | <b><u>Qty</u></b>       |          |
|---|-------------------------|----------|
| <b>Animal Identification: Catch a Predator Kit I (Catalog# 7005)</b>  |                         |          |
| --Sample A, B, C, and Unknown DNA                                     | 60 $\mu$ L for each DNA | included |
| --2X Res. Digest Master Mix (includes restriction enzymes and buffer) | 240 $\mu$ L             | included |
| --DNA Ladder  | 30 $\mu$ L              | included |
| <b>Animal Identification: Catch a Predator Kit II (Catalog# 7006)</b> |                         |          |
| --Sample A, B, C, and Unknown DNA                                     | 60 $\mu$ L for each DNA | included |
| --2X PCR Master Mix (includes primers and enzyme)                     | 270 $\mu$ L             | included |
| --ddH <sub>2</sub> O  | 30 $\mu$ L              | included |
| --DNA Ladder  | 30 $\mu$ L              | included |